

## Product data sheet

### Bellovac ABT accessories

#### 1. Company data

##### 1.1. Manufacturer:

Name Astra Tech AB  
Address Box 14 SE-431 21, Mölndal  
Country Sverige  
Website [www.astratech.com](http://www.astratech.com)

##### 1.2. Distributor:

Name Astra Tech A/S  
Address Husby Allé 19  
Phone no. +45 4374 3377  
Fax no. +45 4371 7865  
E-mail [info.dk@astratech.com](mailto:info.dk@astratech.com)  
Website [www.astratech.dk](http://www.astratech.dk)

##### Quality management systems:

Manufacturer: ISO 9001 / EN 46000

Distributor ISO 9001 / EN 46000

#### 2. Product

2.1 Trade name Bellovac ABT (Transfusion set and Extension line)

##### 2.2. List of product range:

See "Product range" (if Excel file)

"Click here" (if web based)

##### 2.3. More information about the product

See "More information" (if Excel file)

"Click here" (if web based)

#### 3. Product components

(For sterile products including the item package/sterile barrier)

Cardboard:  NO Paper: Laminated virgin paper Glass:  NO

Metal:  NO Latex:  NO Cotton:  NO

Other: \_\_\_\_\_ Organic:  NO

##### 3.2. Plastics:

Types of plastic: Polyamide/Polyethylene, PVC:  YES

If PVC, state plasticiser: DEHP % of PVC: 25-30

##### 3.3. Heavy metals, chlorine, etc.:

Heavy metals:  NO Types: \_\_\_\_\_

Chlorine bleached material:  NO Recycled materials:  NO

Are substances from the Environmental Agency's list of undesirable substances included? : DEHP, cyclohexanone  
State substances: Note: By agreement with Sinerfa, this point takes effect only from 1 June 2003!!

#### 4. Packaging components

(For sterile products **excluding** the item package/sterile barrier)

Cardboard:

Paper:

Glass:

Metal:

Latex:

Cotton:

Organic:

Other: \_\_\_\_\_

#### 4.2. Plastics:

Types of plastic: \_\_\_\_\_

PVC:

If PVC, state plasticiser: \_\_\_\_\_

% of PVC: \_\_\_\_\_

#### 4.3. Heavy metals, chlorine, etc.:

Heavy metals:

Types: \_\_\_\_\_

Chlorine bleached material:

Recycled materials:

State any substances from the Environmental Agency's list of undesirable substances: \_\_\_\_\_

*Note: By agreement with Sinerfa this point takes effect only from 1 June 2003!!*

#### 5. CE marking

5.1. CE mark: No:

Yes:

Class: II a

#### 6. Environmental information

6.1. Does a recognised env. label exist for the product area?

Yes:

No:

6.2. Does the product carry an environmental label?

No:

Yes:

Label: \_\_\_\_\_

6.3. Environmental management systems: (State recognised environmental management systems, if any)

Manufacturer: ISO 14001

Distributor: \_\_\_\_\_

#### 7. Shelf life

Shelf life after date of manufacture / date of sterilisation:

#### 8. Recommended storage and handling

State recommendations Dry environment, normal inhouse temperature

#### 9. Recommendations for waste handling

Recommendation = X	Incineration	Landfill	Recycling	Recycling label
Product	X			
Item package	X			
Inner package			X	X
Transport packaging			X	X

## Bellovac ABT Accessories

### List of product range

Manufacturer product no.	Distributor product no.	EAN no. sales unit	Description	Size	Packaging unit	No. of packs per trsp carton	No. of cartons per pallet	Sterile
REF 68380	REF 68380	1739253200712	Extension line (Jehovah's Witnesses)	All	1			*
REF 68380	REF 68380	5739253200713	Extension line (Jehovah's Witnesses)	All	5	4	20	*
REF 68381	REF 68381	1739253200714	Transfusion set	All	1			*
REF 68381	REF 68381	20739253200715	Transfusion set	All	20	4	20	*

Note: EAN first line units and second line ward box.

#### Instructions:

There must be an EAN no. for each possible sales unit

The description may be the text of the distributor's invoice / delivery note

The size concept varies for different product areas, for instance S-M-L for absorbent pads, measurements, weight and ch for catheters, etc.

A packaging unit is for instance the number of items, rolls, cans, ml or similar.

Indicate sterile products by means of an \*