

Product data sheet

Bellovac

1. Company data

1.1. Manufacturer:

Name Astra Tech AB
 Address Box 14 SE-431 21 Mölndal
 Country Sverige
 Website www.astratech.com

1.2. Distributor:

Name Astra Tech A/S
 Address Husby Allé 19
 Phone no. + 45 4374 3377
 Fax no. + 45 4371 7865
 E-mail info.dk@astratech.com
 Website www.astratech.dk

Quality management systems:

Manufacturer: ISO 9001, EN 46000

Distributor: ISO 9001, EN 46000

2. Product

2.1 Trade name Bellovac

2.2. List of product range:

[See "Product range" \(if Excel file\)](#)

["Click here" \(if web based\)](#)

2.3. More information about the product

[See "More information" \(if Excel file\)](#)

["Click here" \(if web based\)](#)

3. Product components

(For sterile products including the item package/sterile barrier)

Cardboard: NO Paper: Laminated virgin paper Glass: NO

Metal: NO Latex: NO Cotton: NO

Other: Silicone, bariumsulphate, lacquer Organic: NO

3.2. Plastics:

Types of plastic: Polyethylene, Polyamid, ABS, EBA, PP, EVA PVC: YES

If PVC, state plasticiser: DEHP % of PVC: 25-30

3.3. Heavy metals, chlorine, etc.:

Heavy metals: NO Types: _____

Chlorine bleached material: NO Recycled materials: NO

Are substances from the Environmental Agency's list of undesirable substances included? : DEHP

State substances: Note: By agreement with Sinerfa, this point takes effect only from 1 June 2003!!

4. Packaging components

(For sterile products **excluding** the item package/sterile barrier)

Cardboard:

Paper:

Glass:

Metal:

Latex:

Cotton:

Organic:

Other: _____

4.2. Plastics:

Types of plastic: _____

PVC:

If PVC, state plasticiser: _____

% of PVC: _____

4.3. Heavy metals, chlorine, etc.:

Heavy metals:

Types: _____

Chlorine bleached material:

Recycled materials:

State any substances from the Environmental Agency's list of undesirable substances: _____

Note: By agreement with Sinerfa this point takes effect only from 1 June 2003!!

5. CE marking

5.1. CE mark: No:

Yes:

Class: II a

6. Environmental information

6.1. Does a recognised env. label exist for the product area? Yes: No:

6.2. Does the product carry an environmental label? No: Yes: Label: _____

6.3. Environmental management systems: (State recognised environmental management systems, if any)

Manufacturer: ISO 14001

Distributor: _____

7. Shelf life

Shelf life after date of manufacture / date of sterilisation:

8. Recommended storage and handling

State recommendations Dry environment, normal inhouse temperature

9. Recommendations for waste handling

Recommendation = X	Incineration	Landfill	Recycling	Recycling label
Product	X			
Item package	X			
Inner package			X	X
Transport packaging			X	X

Bellovac

List of product range

Manufacturer product no.	Distributor product no.	EAN no. sales unit	Description	Size	Packaging unit	No. of packs per trsp carton	No. of cartons per pallet	Ste- rile
68338	68338	1 7392532 00642 5	Wound drainage system	CH 10	1			*
68338	68338	15 7392532 00643 2	Wound drainage system	CH 10	15	2	20	*
68339	68339	1 7392532 00644 9	Wound drainage system	CH 14	1			*
68339	68339	15 7392532 00645 6	Wound drainage system	CH 14	15	2	20	*
68340	68340	1 7392532 00646 3	Wound drainage system	CH 18	1			*
68340	68340	15 7392532 00647 0	Wound drainage system	CH 18	15	2	20	*

Note: EAN unit (first line) EAN ward box (second line)

Instructions:

There must be an EAN no. for each possible sales unit

The description may be the text of the distributor's invoice / delivery note

The size concept varies for different product areas, for instance S-M-L for absorbent pads, measurements, weight and ch for catheters, etc.

A packaging unit is for instance the number of items, rolls, cans, ml or similar.

Indicate sterile products by means of an *