

## Product data sheet

### Exudrain with trocar

#### 1. Company data

##### 1.1. Manufacturer:

Name Astra Tech AB  
Address Box 14 SE-431 21 Mölndal  
Country Sverige  
Website www.astratech.com

##### 1.2. Distributor:

Name Astra Tech A/S  
Address Husby Allé 19  
Phone no. + 45 43 3377  
Fax no. + 45 43 7865  
E-mail info.dk@astratech.com  
Website www.astratech.dk

##### Quality management systems:

Manufacturer: ISO 9001, EN 46000

Distributor: ISO 9001, EN 46000

#### 2. Product

2.1 Trade name Exudrain

##### 2.2. List of product range:

See "Product range" (if Excel file)

"Click here" (if web based)

##### 2.3. More information about the product

See "More information" (if Excel file)

"Click here" (if web based)

#### 3. Product components

(For sterile products including the item package/sterile barrier)

Cardboard:  Paper: Laminated virgin pa Glass:

Metal: Stainless steel Latex:  Cotton:

Other: Silicone, Bariumsulphate, lacquer Organic:

##### 3.2. Plastics:

Types of plastic: Polyethylene, Polyamid, Polyurethane, ABS, PP, EVA PVC:

If PVC, state plasticiser: DEHP % of PVC: 25-30

##### 3.3. Heavy metals, chlorine, etc.:

Heavy metals:  Types: \_\_\_\_\_

Chlorine bleached material:  Recycled materials:

Are substances from the Environmental Agency's list of undesirable substances included? : DEHP, Cyclohexanone  
State substances: Note: By agreement with Sinerfa, this point takes effect only from 1 June 2003!!

#### 4. Packaging components

(For sterile products **excluding** the item package/sterile barrier)

Cardboard:  Paper:  Glass:   
Metal:  Latex:  Cotton:   
Organic:   
Other: \_\_\_\_\_

#### 4.2. Plastics:

Types of plastic: \_\_\_\_\_ PVC:   
If PVC, state plasticiser: \_\_\_\_\_ % of PVC: \_\_\_\_\_

#### 4.3. Heavy metals, chlorine, etc.:

Heavy metals:  Types: \_\_\_\_\_  
Chlorine bleached material:  Recycled materials:

State any substances from the Environmental Agency's list of undesirable substances: \_\_\_\_\_  
*Note: By agreement with Sinerfa this point takes effect only from 1 June 2003!!*

#### 5. CE marking

5.1. CE mark: No:  Yes:  Class: I I a

#### 6. Environmental information

6.1. Does a recognised env. label exist for the product area? Yes:  No:

6.2. Does the product carry an environmental label? No:  Yes:  Label: \_\_\_\_\_

6.3. Environmental management systems: (State recognised environmental management systems, if any)

Manufacturer: ISO 14001 Distributor: \_\_\_\_\_

#### 7. Shelf life

Shelf life after date of manufacture / date of sterilisation:

#### 8. Recommended storage and handling

State recommendations Dry environment, normal inhouse temperature

#### 9. Recommendations for waste handling

| Recommendation = X  | Incineration | Landfill | Recycling | Recycling label |
|---------------------|--------------|----------|-----------|-----------------|
| Product             | X            |          |           |                 |
| Item package        | X            |          |           |                 |
| Inner package       |              |          | X         | X               |
| Transport packaging |              |          | X         | X               |

## Exudrain with trocar

### List of product range

| Manufacturer product no. | Distributor product no. | EAN no. sales unit | Description                       | Size  | Packaging unit | No. of packs per trsp carton | No. of cartons per pallet | Sterile |
|--------------------------|-------------------------|--------------------|-----------------------------------|-------|----------------|------------------------------|---------------------------|---------|
| 68348                    | 68348                   | 1 7392532 00662 3  | Wound drainage system with trocar | CH 10 | 1              |                              |                           | *       |
| 68348                    | 68348                   | 15 7392532 00663 0 | Wound drainage system with trocar | CH 10 | 15             | 3                            | 20                        | *       |
| 68349                    | 68349                   | 1 7392532 00664 7  | Wound drainage system with trocar | CH 14 | 1              |                              |                           | *       |
| 68349                    | 68349                   | 15 7392532 00665 4 | Wound drainage system with trocar | CH 14 | 15             | 3                            | 20                        | *       |
| 68350                    | 68350                   | 1 7392532 00666 1  | Wound drainage system with trocar | CH 18 | 1              |                              |                           | *       |
| 68350                    | 68350                   | 15 7392532 00667 8 | Wound drainage system with trocar | CH 18 | 15             | 3                            | 20                        | *       |
|                          |                         |                    |                                   |       |                |                              |                           |         |
|                          |                         |                    |                                   |       |                |                              |                           |         |

Note: EAN unit (first line) EAN ward box (second line)

#### Instructions:

There must be an EAN no. for each possible sales unit

The description may be the text of the distributor's invoice / delivery note

The size concept varies for different product areas, for instance S-M-L for absorbent pads, measurements, weight and ch for catheters, etc.

A packaging unit is for instance the number of items, rolls, cans, ml or similar.

Indicate sterile products by means of an \*