

## Product data sheet

### MEDENA Airway suction

#### 1. Company data

##### 1.1. Manufacturer:

Name Astra Tech AB  
Address Box 14 SE-431 21 Mölndal  
Country Sverige  
Website [www.astratech.com](http://www.astratech.com)

##### 1.2. Distributor:

Name Astra Tech A/S  
Address Husby Allé 19  
Phone no. +45 4374 3377  
Fax no. +45 4371 7865  
E-mail [info.dk@astratech.com](mailto:info.dk@astratech.com)  
Website [www.astratech.com](http://www.astratech.com)

##### Quality management systems:

Manufacturer: ISI9001/EN46000

Distributor: ISO 9001/EN46000

#### 2. Product

2.1 Trade name MEDENA-Airway suction

##### 2.2. List of product range:

See "Product range" (if Excel file)

"Click here" (if web based)

##### 2.3. More information about the product

See "More information" (if Excel file)

"Click here" (if web based)

#### 3. Product components

(For sterile products including the item package/sterile barrier)

Cardboard:  NO Paper:  laquered paper Glass:  NO

Metal:  NO Latex:  NO Cotton:  NO

Other: \_\_\_\_\_ Organic:  NO

##### 3.2. Plastics:

Types of plastic: Polypropene/polyethene PVC:  YES

If PVC, state plasticiser: DEHP % of PVC: 30

##### 3.3. Heavy metals, chlorine, etc.:

Heavy metals:  NO Types: \_\_\_\_\_

Chlorine bleached material:  NO Recycled materials:  NO

Are substances from the Environmental Agency's list of undesirable substances included? :  DEHP, Cyclohexanon  
State substances: Note: By agreement with Sinerfa, this point takes effect only from 1 June 2003!!

#### 4. Packaging components

(For sterile products **excluding** the item package/sterile barrier)

Cardboard:

Paper:

Glass:

Metal:

Latex:

Cotton:

Organic:

Other: \_\_\_\_\_

#### 4.2. Plastics:

Types of plastic: \_\_\_\_\_

PVC:

If PVC, state plasticiser: \_\_\_\_\_

% of PVC: \_\_\_\_\_

#### 4.3. Heavy metals, chlorine, etc.:

Heavy metals:

Types: \_\_\_\_\_

Chlorine bleached material:

Recycled materials:

State any substances from the Environmental Agency's list of undesirable substances: \_\_\_\_\_

*Note: By agreement with Sinerfa this point takes effect only from 1 June 2003!!*

#### 5. CE marking

5.1. CE mark: No:

Yes:

Class: I

#### 6. Environmental information

6.1. Does a recognised env. label exist for the product area? Yes:  No:

6.2. Does the product carry an environmental label? No:  Yes:  Label: \_\_\_\_\_

6.3. Environmental management systems: (State recognised environmental management systems, if any)

Manufacturer: ISO14001

Distributor: \_\_\_\_\_

#### 7. Shelf life

Shelf life after date of manufacture / date of sterilisation:

#### 8. Recommended storage and handling

State recommendations Dry environment, normal inhouse temperature

#### 9. Recommendations for waste handling

Recommendation = X	Incineration	Landfill	Recycling	Recycling label
Product	X			
Item package	X			
Inner package			X	X
Transport packaging			X	X

## MEDENA Airway suction

### List of product range

Manufacturer product no.	Distributor product no.	EAN no. sales unit	Description	Size	Packaging unit	No. of packs per trsp carton	No. of cartons per pallet	Ste- rile
REF 68106	REF 68106	1739253200490	Airway suction,Coupling nipple/cone	Length 2m	1			*
REF 68106	REF 68106	50739253200491	Airway suction,Coupling nipple/cone	Length 2m	50	2	20	*
REF 68114	REF 68114	1739253200506	Airway suction, Coupling nipple/cone	Length 3,5 m	1			*
REF 68114	REF 68114	35739253200507	Airway suction, Coupling nipple/cone	Length 3,5 m	35	2	20	*
REF 68207	REF68207	1739253200596	Airw.suct.,Coupl.nipple/cone,fing. cont.h	Length 2 m	1			*
REF 68207	REF68207	50739253200597	Airw.suct.,Coupl.nipple/cone,fing. cont.h	Length 2 m	50	2	20	*

Note: EAN unit (13 digits) EAN Transport box (14 digits)

Instructions:

There must be an EAN no. for each possible sales unit

The description may be the text of the distributor's invoice / delivery note

The size concept varies for different product areas, for instance S-M-L for absorbent pads, measurements, weight and ch for catheters, etc.

A packaging unit is for instance the number of items, rolls, cans, ml or similar.

Indicate sterile products by means of an \*